

## Dress Sling Padded Roll

## ST2A

- Easy to apply padded sling
- Underarm padded roll supplied for comfort
- Built in waist support with clunk-click fastening
- Good access for hygiene purposes
- Good upper body tone required
- Full risk assessment required for suitability
- Colour coded tapes
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:



100% Polyester

		A	B	C	D
XXXS	ST2AXXS	38 - 58cm	18cm or over	53cm	17½cm
XXS	ST2AXSM	59 - 77cm	20cm or over	60cm	19cm
XS	ST2AXSM	77 - 87cm	22cm or over	73cm	21cm
S	ST2ASM	87 - 107cm	25cm or over	83cm	25cm
M	ST2AMED	97 - 117cm	27cm or over	93cm	27cm
L	ST2ALRG	112 - 132cm	28cm or over	102cm	29½cm
XL	ST2AXLRG	125 - 145cm	30cm or over	112cm	33cm
XXL	ST2AXXLRG	138 - 158cm	34cm or over	122cm	38cm

