



THERMOGUN® TG2

Medical infrared forehead thermometer with digital display

- Hygienic non-contact temperature reading
- Medical grade sensor
- Fast and accurate with low tolerance range
- Simple and easy to use with bright LCD display – ideal for night time
- 3 different display backlights for individual temperature ranges: Green, Orange, Red
- Alarm function
- Memorization of the last 32 measurements
- Tested and approved by TÜV

⊕ NON CONTACT ⊕ ACCURATE ⊕ INTELLIGENT



THERMOGUN TG2 - MEDICAL INFRARED THERMOMETER

Art. Nr. 64001

EAN 4260041686045

FEATURES & SPECIFICATIONS

- Hygienic non-contact temperature reading
- Medical grade sensor
- Suitable for medical use
- Fast and accurate with low tolerance range
- Simple and easy to use with bright LCD display – ideal for night time
- Three different display backlights for individual temperature ranges: Green, Orange, Red
- Memorization of the last 32 measurements
- Alarm function
- Auto standby function
- Measuring distance: < 2.5cm
- Detection range
Body: 32.0°C - 42.9°C (Accuracy ±0.2°C)
Objects: 0°C - 100.0°C (Accuracy ±1.0°C)
- Operating environment
Temperature: 10°C - 40°C
Relative humidity: ≤ 85%
Atmospheric pressure: 70kPa - 106kPa
- Storage conditions
Temperature: -20°C - +55°C
Relative humidity: ≤ 93%
Atmospheric pressure: 70kPa - 106kPa
- Power supply: DC3V - 2xAA batteries (not included)
- Size: 160 x 95 x 46 mm
- Weight: 144g

Suitable for medical use:
at home, in public places,
medical practices or hospitals



Technical changes and mistakes reserved.
The product images and colors may differ slightly.



With 3 different display backlights for individual temperature ranges:

